THUASNE

Ligaflex® Immo

fr	Attelle d'immobilisation du poignet évolutive	5
en	Modular wrist immobilisation splint	6
de	Modulare Handgelenksorthese	8
nl	Brace voor immobilisatie en progressieve opbouw van de pols	10
it	Tutore di immobilizzazione evolutiva per polso	12
es	Órtesis evolutiva de inmovilización de la muñeca	14
pt	Tala de imobilização do punho, evolutiva	16
da	Dynamisk skinne til immobilisering af håndled	18
fi	Modulaarinen ranteen immobilisaatiotuki	20
sv	Anpassningsbar handledsskena för immobilisering	21
el	Προσαρμόσιμος várθηκαç ακινητοποίησηç καρπού	23
cs	Modulární imobilizační ortéza zápěstí	25
pl	Orteza nadgarstkowo-ręczna, unieruchamiająca z możliwością adaptacji	27
lv	Regulējamā šina plaukstas locītavas imobilizēšanai	29
lt	Pritaikomas riešo imobilizavimo įtvaras	31
et	Reguleeritav randmelahas	33
sl	Prilagodljiva opornica za imobilizacijo zapestja	34
sk	Modulárna imobilizačná ortéza zápästia	36
hu	Moduláris csuklórögzítő sín	38
bg	Шина за имобилизация на китка, адаптируема	40
ro	Atelă ajustabilă de imobilizare a încheieturii mâinii	42
ru	Шина для имобилизации ревматоидного запястья	44
hr	Udlaga za imobilizaciju ručnog zgloba, produljiva	46
zh	有适应能力的手腕固定夹板	47
ar	جبيرة لتثبيت الرسغ	53

	1	13 - 15 cm
	2	15,5 - 17,5 cm
	3	18 - 20 cm
	4	20,5 - 23 cm

fr	Tour de poignet
en	Wrist circumference
de	Handgelenkumfang
nl	Polsomtrek
it	Circonferenza del polso
es	Contorno de la muñeca
pt	Perímetro do punho
da	Håndleddets omkreds
fi	Ranteen ympärysmitta
sv	Handledens omkrets
el	Περίμετρος καρπού
cs	Obvod zápěstí
pl	Obwód nadgarstka
lv	Plaukstas locītavas apkārtmērs
it	Riešo apmītis
et	Randme ümbermõõt
sl	Obseg zapestja
sk	Obvod zápästia
hu	Csukló kerülete
bg	Обиколка на китката
ro	Circumferința încheieturii
ru	Обхват запястья
hr	Obujam ručnog zgloba
zh	腕围
ar	محيط الرسغ

fr	Immobilisation
en	Immobilisation
de	Ruhigstellung
nl	Immobilisatie
it	Immobilizzazione
es	Inmovilización
pt	Imobilização
da	Støtte
fi	Immobilisoi
sv	Immobilisering
el	Ακινητοποίηση
cs	Znehybnění
pl	Unieruchomienie
lv	Imobilizācija
lt	Imobilizavimas
et	Liikumatus muutmine
sl	Imobilizacija
sk	Imobilizácia
hu	Rögzítés
bg	Обездвижване
ro	Imobilizare
ru	Иммобилизация
hr	Imobilizacija
zh	固定
ar	التثبيت

MD

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MODULAR WRIST IMMOBILISATION SPLINT

Description/Destination

This device is a wrist splint providing joint immobilisation.

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

polyester - cotton - elastane - polyurethane - polyoxymethylene - aluminium - rubber.

Propriétés/Mode of action

Immobilisation of the wrist via a rigid, malleable and anatomical palmar stay and a dorsal stay.

Aerated lining.

Possibility of immobilising the wrist in 3 different positions.

Ligaflex[®] Immo is supplied with a pair of reinforcements / stays enabling immobilisation of the wrist in the extended position (radiocarpal extension).

It is possible to order a kit containing 2 pairs of reinforcements separately;
- Reinforcements immobilizing the wrist in the flexed position (radio carpal flexion).
- Reinforcements immobilizing the wrist in the flexed and inclined position (palmar flexion and cubital inclination).

Indications

Wrist sprain.

Non surgical wrist fracture.

Non surgical scaphoid fracture.

Post-operative or post-traumatic wrist immobilisation.

Wrist tendinopathies.

Lateral epicondylitis.

Rheumatoid arthritis.

Conservative treatment of carpal tunnel syndrome.

Neuro-muscular impairments (contracture, hand drop or spasticity...).

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

The choice of reinforcements should be determined on the basis of the prescription of the health professional supplying the device. Do not change the reinforcements without first consulting your professional.

In the event of a palmar injury, do not use the device in the extended wrist position.

In the event of an extension injury, do not use the flexed wrist positions.

Extended wrist position (radiocarpal extension).

Flexed wrist position (radiocarpal flexion).

Flexed wrist position (radiocarpal flexion) and cubital inclination.

To insert the reinforcements:

Unfasten the self-fastening part at the base of the splint ❶. Remove the reinforcements.

Select the pair of reinforcements on the basis of the condition treated.

Slip the stay inside the reinforcement cover, making sure that the pictogram is visible. ❷

These pictograms refer to palmar stays (under the hand).

These pictograms refer to dorsal stays (over the hand).

Open the splint and position the wrist inside the product.

Fasten the straps, following the order in which the straps are numbered (1 point, 2 points, 3 points).

Care

Product can be washed in accordance with the instructions shown on this leaflet and on the label. Remove metal stays before washing. Replace in the original location before next use. Close the self-fastening tabs before washing. Machine washable at 40°C (delicate programme). Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not dry clean. Do not tumble-dry.

Do not iron. Squeeze out excess water. Dry flat. Dry away from any direct heat source (radiator, sun, etc.).

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2027601 (2021-04)